


CLD DOGBONE



TyeTec® CLD Dogbones wurden speziell für die Verwendung mit TyeTec® padeyes Typ „DL“ und „DF“ entworfen. Ihre Form ist dem Innenradius des Padeyes angepasst um die Kraft sicher einzuleiten. Die Einschnürung in der Mitte hält den Dogbone sicher an seinem Platz. Für manche Anwendungen kann es sinnvoll sein, den Dogbone mit ein paar Stichen und einem Takling zu sichern.

TyeTec® dogbones were specially designed for the use with TyeTec® padeyes. Their shape is adjusted to the inner radius of the padeye to transfer power safely. The construction in the center keeps the dogbone in place. For some applications, it may be useful to secure the dogbone with a few stitches and a takling.

Product nr.	L	D	weight	applications
CLD06	12	8	3	DL06...
CLD08	20	10	12	DL08..., DF 08..., DM08ALS, SP 08, W45 single Loop
CLD10	25	10	25	DL10..., DF10..., DM10ALS, SP10, W60 single Lopp
CLD12	30	15	40	DL12...,DF12..., DM12ALS, SP12, W80 single Loop
CLD16	40	20	90	DL16..., DM16ALS, W 100 single Loop
CLD24	60	28	250	DL24..., W120 single Loop